

# Start M3 Hemi

## Ideal for hemiplegia fittings



The Start M3 Hemi is designed for hemiplegia fittings.

Thanks to its small anterior seat height and a wider foot space, the user can propel the wheelchair with his or her sound foot, provided that the foot rest has been removed beforehand.

A comprehensive selection of options is available for the Start M3 Hemi, such as the one-arm drive, knee lever wheel lock for one-handed operation and channel forearm pad for the affected arm, rounding out the fitting.



- The hemi kit serves to widen the foot space and permits a smaller anterior seat height.



- The one-arm drive with telescoping axle allows independent mobility.



- The channel forearm pad positions and stabilises the affected arm.

# Options

---



• Rotation adapter for channel forearm pads



• Silicone push ring cover offers enhanced grip



• Double push rings permit one-hand operation and ensure directional stability



• Functional crutch holder is ideal for all walking aids

# Product Features

- Wider foot space for small anterior seat height
- Wheelbase extension adjustment possible
- Aluminium tube frame
- Modular design
- All main component groups are attached with screws
- Durable upholstery
- Narrow folding size with folding fixture
- Sturdy one-and-half fold crossbrace
- Optimised knee lever wheel lock



## Basic Technical Data

Load capacity	125 kg
Weight	From 15.5 kg
Seat width	38 cm - 50.5 cm
Seat depth	40 cm, 43 cm, 46 cm
Back height	40 cm - 47.5 cm
Seat height	Front 38 cm - 51 cm Rear 38 cm - 49 cm
Lower leg length	27 cm - 47 cm