# VICAIR® CENTRE RELIEF

# **OFFLOADING PRESSURE IN THE HIGH RISK AREA**



# For completely offloading of the rectal area

# **Adjustable** non-contact area

## Stable positioning

The Vicair Centre Relief is the best choice if you need to completely offload the rectum perineum and coccyx area. By simply removing SmartCells from the centre circle compartment you create a non-contact area. This important feature makes a great difference for persons with (for instance) localized deep pressure ulcers, rectum or tailbone problems. Or after perineal surgery.

The Vicair Centre Relief helps you to sit comfortably and safe, even in the most painful situations.

#### PRODUCT PERFORMANCE Skin protection / Anti decubitus Pressure redistribution .... Shear force reduction

Microclimate control	
Impact & vibration damping	

### Positioning

i usitioning	
Lateral stability	
Frontal stability	
Sliding reduction	
Allows manual adjustments	
Allows easy transfers	

Maintenance & Safety	
Durability	
Reliability	
Functionality when misused	
Low maintenance level	
Level of hygiene	
,0	

Functionality & Application	
Comfort	
Low cushion weight	
Available size range	
Maximum user weight	
Low seat to floor height	
Auto adjustments	




## **USED WITH**

Rectum problems, tailbone problems, localized deep pressure ulcers, after perineal surgery or difficult childbirth

TECHNICAL SPECIFICATIONS		
Height	10 cm	
Average weight	750 grams	
Max. user weight	250 kg	

## **CUSHION SIZES** Standard On request



070069-ACR-EN-REV04/2019