

SIZE SELECTION CHART FOR YOKE SLINGS

This document is intended as a guide only. Correct selection for an individual patient remains the responsibility of the healthcare professional or case manager.

It is critical to select the correct size of sling to match the individual patient. Typically there is a degree of overlap between each


sling size within a particular model. If the patient falls between two sizes, it is generally better to select the smaller size of sling, unless you are aware the patient is likely to grow.

Please note, however, that selecting a sling involves more than simply measuring a patient's size.

To be confident a sling meets individual requirements and is the most appropriate sling available, the following basic factors should be taken into consideration:

- patient information
- lifting task and environment
- carer information
- other equipment information

For more details, please refer to our ProSling Range Guide available to download from www.novis.com.au



Yoke General Purpose	S Shoulder width mm																	C Coccyx to Shoulders mm																																							
	500-520	520-540	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080	1080-1100	1100-1120	1120-1140	1140-1160	1160-1180	1180-1200	1200-1220	1220-1240	1240-1260	1260-1280	1280-1300	440-460	460-480	480-500	500-520	520-540	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780
SL147181	Small															Small																																									
SL147281	Medium																	Medium																																							
SL147381	Large																					Large																																			

Yoke General Purpose with Head Support	S Shoulder width mm																	H Coccyx to top of head mm (Sitting Height)																																																			
	500-520	520-540	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080	1080-1100	1100-1120	1120-1140	1140-1160	1160-1180	1180-1200	1200-1220	1220-1240	1240-1260	1260-1280	1280-1300	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080	1080-1100	1100-1120	1120-1140	1140-1160	1160-1180	1180-1200	1220-1240
SL146181	Small																	Small																																																			
SL146281	Medium																	Medium																																																			
SL146381	Large																	Large																																																			
SL146481	Extra Large																	Extra Large																																																			
SL146581	Extra Extra Large													Extra Extra Large																																																							

Yoke Hygiene with Head Support	S Shoulder width mm																	T Sacrum to top of head mm (Sitting Height)																																																			
	500-520	520-540	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080	1080-1100	1100-1120	1120-1140	1140-1160	1160-1180	1180-1200	1200-1220	1220-1240	1240-1260	1260-1280	1280-1300	500-520	520-540	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080
SL246181	Small																	Small																																																			
SL246281	Medium																	Medium																																																			
SL246381	Large																	Large																																																			
SL246481	Extra Large																	Extra Large																																																			
SL246581	Extra Extra Large													Extra Extra Large																																																							

Yoke Hygiene	S Shoulder width mm																	W Waist circumference mm																																																				
	500-520	520-540	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080	1080-1100	1100-1120	1120-1140	1140-1160	1160-1180	1180-1200	1200-1220	1220-1240	1240-1260	1260-1280	1280-1300	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080	1080-1100	1100-1120	1120-1140
SL267181	Small																	Small																																																				
SL267281	Medium																	Medium																																																				
SL267381	Large																	Large																																																				

Disposable General Purpose with Head Support	S Shoulder width mm																	H Coccyx to top of head mm (Sitting Height)																
	500-520	520-540	540-560	560-580	580-600	600-620	620-640	640-660	660-680	680-700	700-720	720-740	740-760	760-780	780-800	800-820	820-840	840-860	860-880	880-900	900-920	920-940	940-960	960-980	980-1000	1000-1020	1020-1040	1040-1060	1060-1080	1080-1100	1100-1120	1120-1140		
SL156281	Medium																	Medium																
SL156381	Large																	Large																



SIZE SELECTION CHART FOR PIVOT SLINGS

Pivot Hygiene	S Shoulder width mm																																						
	500 - 520	520 - 540	540 - 560	560 - 580	580 - 600	600 - 620	620 - 640	640 - 660	660 - 680	680 - 700	700 - 720	720 - 740	740 - 760	760 - 780	780 - 800	800 - 820	820 - 840	840 - 860	860 - 880	880 - 900	900 - 920	920 - 940	940 - 960	960 - 980	980 - 1000	1000 - 1020	1020 - 1040	1040 - 1060	1060 - 1080	1080 - 1100	1100 - 1120	1120 - 1140	1140 - 1160	1160 - 1180	1180 - 1200	1200 - 1220	1220 - 1240	1240 - 1260	1260 - 1280
SL247191	Small																																						
SL247291					Medium																																		
SL247391					Large																																		

Pivot General Purpose with Head Support	S Shoulder width mm																																						
	500 - 520	520 - 540	540 - 560	560 - 580	580 - 600	600 - 620	620 - 640	640 - 660	660 - 680	680 - 700	700 - 720	720 - 740	740 - 760	760 - 780	780 - 800	800 - 820	820 - 840	840 - 860	860 - 880	880 - 900	900 - 920	920 - 940	940 - 960	960 - 980	980 - 1000	1000 - 1020	1020 - 1040	1040 - 1060	1060 - 1080	1080 - 1100	1100 - 1120	1120 - 1140	1140 - 1160	1160 - 1180	1180 - 1200	1200 - 1220	1220 - 1240	1240 - 1260	1260 - 1280
SL146191					Small																																		
SL146291									Medium																														
SL146391													Large																										

This document is intended as a guide only. Correct selection for an individual patient remains the responsibility of the healthcare professional or case manager.

It is critical to select the correct size of sling to match the individual patient. Typically there is a degree of overlap between each


sling size within a particular model. If the patient falls between two sizes, it is generally better to select the smaller size of sling, unless you are aware the patient is likely to grow.

Please note, however, that selecting a sling involves more than simply measuring a patient's size.

To be confident a sling meets individual requirements and is the most appropriate sling available, the following basic factors should be taken into consideration:

- patient information
- lifting task and environment
- carer information
- other equipment information

For more details, please refer to our ProSling Range Guide available to download from www.novis.com.au



H Coccyx to top of head mm (Sitting Height)																																								
	540 - 560	560 - 580	580 - 600	600 - 620	620 - 640	640 - 660	660 - 680	680 - 700	700 - 720	720 - 740	740 - 760	760 - 780	780 - 800	800 - 820	820 - 840	840 - 860	860 - 880	880 - 900	900 - 920	920 - 940	940 - 960	960 - 980	980 - 1000	1000 - 1020	1020 - 1040	1040 - 1060	1060 - 1080	1080 - 1100	1100 - 1120	1120 - 1140										
SL146191					Small																																			
SL146291													Medium																											
SL146391																	Large																							



SIZE SELECTION CHART FOR GAIT TRAINING VEST

Gait Training Vest	W Waist circumference mm																																															
	520 - 540	540 - 560	560 - 580	580 - 600	600 - 620	620 - 640	640 - 660	660 - 680	680 - 700	700 - 720	720 - 740	740 - 760	760 - 780	780 - 800	800 - 820	820 - 840	840 - 860	860 - 880	880 - 900	900 - 920	920 - 940	940 - 960	960 - 980	980 - 1000	1000 - 1020	1020 - 1040	1040 - 1060	1060 - 1080	1080 - 1100	1100 - 1120	1120 - 1140	1140 - 1160	1160 - 1180	1180 - 1200	1200 - 1220	1220 - 1240	1240 - 1260	1260 - 1280	1280 - 1300	1300 - 1320								
SL346191													Small																																			
SL346291																					Medium																											
SL346391																																	Large															

This document is intended as a guide only. Correct selection for an individual patient remains the responsibility of the healthcare professional or case manager.

It is critical to select the correct size of sling to match the individual patient. Typically there is a degree of overlap between each

sling size within a particular model. If the patient falls between two sizes, it is generally better to select the smaller size of sling, unless you are aware the patient is likely to grow.

Please note, however, that selecting a sling involves more than simply measuring a patient's size.

To be confident a sling meets individual requirements and is the most appropriate sling available, the following basic factors should be taken into consideration:

- patient information
- lifting task and environment
- carer information
- other equipment information

For more details, please refer to our ProSling Range Guide available to download from www.novis.com.au

