

The QUADTRO SELECT® and CONTOUR SELECT® cushions feature four independent compartments of air cells for an individualised fit, each controlled by our exclusive ISOFLO® Memory Control.

ROHO's ISOFLO® Memory Control, located at the front of all Select Series cushions, allows air to flow through the cushion in two ways:



Air can flow freely throughout the entire cushion, similar to a ROHO Single Compartment Cushion.



Air becomes locked into four separate quadrants, offering front-to-back and side-to-side positioning capabilities.

Every Select Series cushion sets the standard for skin and soft tissue protection, without compromising stability, positioning, comfort or ease of use.

Featuring DRY FLOATATION® air cells, the Select Series assists in pressure injury prevention and treatment.