## Start M2S - Start 4 standard and Start 4 adjustable model

In order to ensure the rapid availability of our top model Start M2, there are many stock versions, which we refer to as M2S. The "S" stands in this case for "stock". For you this means that these models are available quickly thanks to the optimum logistics chain and can thus be used quickly by patients.





## Technical data

	M2S Standard	M2S Adjustable
User weight	130 kg	130 kg
Seat width :	405 / 455 / 505 mm	380 / 405 / 430 / 455 / 480mm
Seat depth :	405 / 430 / 455 mm	405 / 430 / 455 mm
Back height :	400 / 425 / 450 / 475 mm	400 / 425 / 450 / 475 mm
Anterior seat height	430 - 520 mm	430 - 520 mm
Posterior seat height	410 - 505 mm	410 - 505 mm
Seat angle	0 – 5°	0 – 5°
Lower leg length	280 – 500 mm (continuously adjustable)	280 – 500 mm (continuously adjustable)
Rear wheels	24" standard	24" standard
Caster wheels	8" PU (200x50mm)	6" PU
Side Panel with armrest	Height adjustable (with tools)	Height adjustable (10cm)(push button)
Frame colour	Silver Metallic	Silver Metallic
Back	With Lumbar curve	Straight
Back upholstery	Standard	Adjustable economy
Push handles	Standard	Standard
Seat upholstery	Standard	Standard

## **Product highlights**

- Stock version of the fully configurable Start M2
- Patented, tool-free, quick seat depth adjustment (50 mm)
- Very low weight with full Ottobock quality (from 14.7 kg total weight, from 9.6 kg transport weight\*)
- Antibacterial push handles\*\*
- Excellent component compatibility
- · New adjustable economy back upholstery
- Short delivery times for stock versions (quick availability)
- Aluminium frame
- · Frame colour: silver metallic
- · Side panel with armrest, height can be adjusted (without tools) or set (with tools) and adjustable depths

<sup>\*</sup> The transport weight is the total weight of the wheelchair minus the weight of the rear wheels, footrests and seat cushion.

<sup>\*\*</sup> The antibacterial effect of the push handles has been verified by an independent institute.